

Niagara Falls
Artist | John Henry Twachtman

The logo for UNIS, featuring a stylized white 'U' with a horizontal line through it, followed by the letters 'NIS' in a bold, sans-serif font.

Power Injectable PICC

- Large lumen Diameter
- Reverse Tapered Design
- Easily identification


GENOSS
For Patients & Doctors

www.genoss.com | 문의 031-888-5100



Power Injectable PICC


PERIPHERALLY INSERTED CENTRAL CATHETER



FLOW RATE COMPARISON

Lumen	Catheter Size (Fr)	Gravity Flow Rate(ml/h)				Max Flow Rate(ml/s)	
		Unis		B Co.		Unis	B Co.
Single	5	1,312		1,185		5	5
Dual	5	660		578			
Triple	6	17G	1,369	17G	1,163	6	6
		19G	463	19G	275		


LARGE LUMEN SIZE



ACCESSORY


1. Peel away sheath

- Easy-to-thread funnel system
- Patented locking system



2. Better Radiopaque Guide wire

- 2 type of length (70/120cm)
- Improvement of radiopacity




Gold coated

Marker

3. Better Visible Needle

- Dimpling method
- Improvement of visibility



Genoss Company A Company B Company C